

POWER

WHEELCHAIR



STAR

STEEL

/ TYPE	Electric Wheelchair
/ CONTROLLER	PG/Dynamic/Haozhi
/ MOTOR	2x300/340W, electromagnetic brake
/ MAX SPEED	6km/h(300W), 10km/h(340W)
/ TRAVEL MILEAGE	30km(38Ah), 45km(52Ah)
/ BATTERY	Lead-acid batteries 24V 38/52Ah
/ FRONT WHEEL	8"x2"/9"x3"
/ REAR WHEEL	12.5"x2"/14"x3"
/ TOTAL MASS(BATTERY INCLUDED)	90kg (38Ah), 104.5kg (52Ah)
/ WEIGHT OF BATTERY	20kg (38Ah), 30kg (52Ah)
/ MAXIMUM SAFE SLOPE	10°
/ MINIMUM TURNING RADIUS	900mm
/ MIN. BRAKING DISTANCE AT MAX. SPEED	6km/h: 1000mm, 10km/h: 1700mm
/ MAXIMUM HEIGHT OF OBSTACLE	50 mm
/ SUSPENSION	Rear Suspension



Product Features



Multi-angle adjustment



Rear lights



Dynamic/PG controller

ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Height Footplate (cm)	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
STAR	24V 38Ah	108	58	94	13~20	43~53	43~49	*	45	136	8"/9"	12.5"/14"	*	*
	24V 52Ah	108	58	94	13~20	43~53	43~49	*	45	136	8"/9"	12.5"/14"	*	*

EASY WAY PLUS

ALUMINUM

/ FRAME	Aluminum, powder coating
/ CONTROLLER	PG/China
/ MOTOR	180W×2, brush motor
/ MAX SPEED	7 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lithium battery 24V 20Ah, 30Ah optional
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12" pneumatic/PU, aluminum alloy
/ ARMREST	Flip back armrest, PU armpad
/ FOOTREST	Flip up, with aluminum pad
/ SEAT	Breathable backrest, thick cushion
/ PACKAGE	95pcs/20"GP, 240pcs/40"HQ
/ SPECIAL	Max. 30Ah lithium battery; Electromagnetic brake; Safety belt; Backrest height adjustable; Tool-free detachable battery box



Product Features



Automatic body folding
with one touch



Rear-foldable armrests on both sides
Controller left and right interchangeable



Intelligent electromagnetic brake

ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
EASY-WAY PLUS		105	67	94	39	45	42	43	51	135	8"	12"	40.5/36	72×43×88

PIONEER

STEEL ALUMINUM

/ FRAME	Steel/Aluminum, powder coating
/ CONTROLLER	PG/China
/ MOTOR	200W×2, brush motor
/ MAX SPEED	7 km/ h
/ TRAVEL MILEAGE	20 km
/ BATTERY	Lead acid 2×12V 20Ah, 32Ah/36Ah optional
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12.5" pneumatic/PU, aluminum alloy
/ ARMREST	Flip back slope armrest, PU armpad, height adjustable optional
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Tension nylon with sponge, thick cushion
/ PACKAGE	75pcs/ 20" GP, 185pcs/ 40" HQ
/ SPECIAL	Backrest half folding; Safety belt; Electromagnetic brake



Product Features



Backrest half folding



Safety belt



Optional: mobile phone holder, cup holder

ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
PIONEER	Steel	108	63	97	34	43-48	43	45	52	135	8"	12.5"	65.6/61.6	87×38×81
	Aluminum	108	63	97	34	43-48	43	45	52	135	8"	12.5"	58.0/50.3	87×38×81

RECLINER

ALUMINUM

/ FRAME	Aluminum, powder coating, foldable
/ CONTROLLER	PG/China
/ MOTOR	340W×2, brush motor
/ MAX SPEED	8-9 KM/H
/ TRAVEL MILEAGE	20 KM
/ BATTERY	Lead acid 2×12V 32Ah
/ FRONT WHEEL	9" PU tire
/ REAR WHEEL	14" solid tire
/ ARMREST	Height adjustable armrest, PU armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Tension nylon with sponge, thick cushion
/ PACKAGE	75pcs/20"GP, 190pcs/40"HQ
/ SPECIAL	Headrest; Elevating legrest; Anti-tipper wheels; Safety belt; Reclining backrest(max. 130°); Footplate height and angle adjustable; Armrest height adjustable(150-170 mm);Seat depth and width adjustable; Anti-tipper wheels



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded (cm)	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
RECLINER		107	65	97	36	44-50	47.5-50	44	53	135	9"	14"	71.7/66.7	87×45×83

SMART

ALUMINUM

/ FRAME	Aluminum, powder coating
/ CONTROLLER	China
/ MOTOR	150W×2, brushless motor
/ MAX SPEED	6 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lithium battery, 12Ah
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12" pneumatic/PU, aluminum alloy
/ ARMREST	Height adjustable armrest, PU armpad
/ FOOTREST	Detachable with angle adjustable footplate
/ SEAT	Tension nylon backrest
/ PACKAGE	Breathable cushion seat
/ SPECIAL	Safety belt; Backrest half folding



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded (cm)	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SMART Kid		100	63	95	49	38.5	41	48	50	120	8"	12"	29.4/24.9	85×45×84
SMART		100	65	95	41	45	41	48	50	120	8"	12"	29.4/24.9	85×45×84

ECONOMY

STEEL

/ FRAME	Steel, powder coating
/ CONTROLLER	China
/ MOTOR	250W×2, brush motor
/ MAX SPEED	7 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lead acid 2x12V 12Ah, 20Ah optional
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	16" pneumatic, steel alloy
/ ARMREST	Flip back slope armrest, PVC armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Tension nylon with sponge, breathable cushion
/ PACKAGE	80pcs/20"GP, 210pcs/40"HQ
/ SPECIAL	Safety belt; Hand brakes



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded (cm)	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
ECONOMY		106	64	94	37	43-48	44	41	53	135	8"	16"	51.0/43.0	87x38x99

JUMPER

ALUMINUM

/ FRAME	Aluminum, powder coating
/ CONTROLLER	PG/China
/ MOTOR	340W×2, brush motor
/ MAX SPEED	9 KM/H
/ TRAVEL MILEAGE	20 KM
/ BATTERY	Lead acid 2x12V 32Ah, 20Ah /36Ah optional
/ FRONT WHEEL	9" PU tire, 10" optional
/ REAR WHEEL	14" pneumatic/solid, aluminum alloy
/ ARMREST	Flip back height adjustable slope armrest, PU armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Tension nylon with sponge, thick cushion
/ PACKAGE	75pcs/20"GP, 160pcs/40"HQ
/ SPECIAL	Phone holder; Cup holder; Safety belt; Electromagnetic brake;



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded (cm)	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height (cm)	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
JUMPER		107	65	97	36	43-48	44	44	53	135	9"	14"	81.0/62.0	89x42x104

SPLIT

ALUMINUM

/ FRAME	Aluminum, powder coating
/ CONTROLLER	PG/ China
/ MOTOR	200W×2, brush motor
/ MAX SPEED	7 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lithium 24V 20Ah (Lead acid 20Ah battery optional)
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12.5" pneumatic/PU, aluminum alloy
/ ARMREST	Flip back slope armrest, PU armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Breathable backrest and cushion
/ PACKAGE	95pcs/20"GP, 245pcs/40"HQ
/ OPTION	Height adjustable armrest; Color upholstery
/ SPECIAL	Backrest half folding; Safety belt; Electromagnetic brake



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SPLIT		108	63	97	34	45	44	44	53	135	8"	12.5"	40.5/33.5	87x38x81

SPLIT II

STEEL

/ FRAME	Steel, powder coating
/ CONTROLLER	PG/China
/ MOTOR	250W×2, brush motor
/ MAX SPEED	7 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lead acid 2×12V, 20Ah
/ FRONT WHEEL	8"×2", PU tire
/ REAR WHEEL	12.5" pneumatic/PU, steel spoke
/ ARMREST	Flip back slope armrest, PU armpad
/ FOOTREST	Detachable footrest with heel support
/ SEAT	Tension nylon with sponge, thick cushion
/ SPECIAL	Electromagnetic brake; Safety belt; Backrest half folding



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SPLIT II	15"	108	63	97	46	38	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81
	18"	108	63	97	46	45	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81
	20"	108	63	97	46	50	44	44	53	135	8"	12.5"	37.0/35.0	87x48x81

ENGULM

STEEL

/ FRAME	Steel, powder coating
/ CONTROLLER	PG/Dynamic
/ MOTOR	340W×2, brush motor
/ MAX SPEED	9 KM/H
/ TRAVEL MILEAGE	25 KM
/ BATTERY	Lead acid, 24V 50Ah (36Ah optional)
/ FRONT WHEEL	9", pneumatic/ PU tire
/ REAR WHEEL	14", pneumatic tire, aluminum alloy
/ ARMREST	Height adjustable armrest, wide PU armpad
/ FOOTREST	Detachable footrest with heel support, width adjustable
/ SEAT	Tension nylon backrest, thick cushion
/ PACKAGE	30pcs/20"GP, 66pcs/40"GP
/ SPECIAL	Powered seat tilting and back reclining; Headrest; Electromagnetic brake; Light system; Seat width adjustable; Seat depth adjustable



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
ENGULM		114.3	62.5	99.5	/	43-50	33-50	49.5	50	135	9"	14"	110.0/98.0	102x77x97

SWIFTY

STEEL

/ FRAME	Steel, powder coating
/ CONTROLLER	PG/China
/ MOTOR	240W×2, brush motor
/ MAX SPEED	6 KM/H
/ TRAVEL MILEAGE	15 KM
/ BATTERY	Lead acid 2×12V 12Ah, lithium battery 20Ah optional
/ FRONT WHEEL	6", black PU tire
/ REAR WHEEL	8", black PU tire
/ ARMREST	Flip back armrest, PU armpad
/ FOOTREST	Flip up aluminum footrest, height adjustable
/ BACKREST	Soft PU leather, long-lasting, easy cleaning
/ SEAT	PU leather, forward and backward adjustable
/ PACKAGE	75pcs/20"GP, 190pcs/40"HQ
/ SPECIAL	Seat width adjustable; Detachable into 3 pieces



ITEM NO.		Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to Seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G./N. Weight (kg)	Carton Size (cm)
SWIFTY		90	53	93	/	46	34	40	53	120	6"	8"	52.5/46.5	82x62x62